



M480 Embedding Cassettes

Made of acetal


Disposable plastic tissue cassettes are suitable for holding and identifying tissue samples in processing, embedding, and sectioning procedures. The cassettes fit securely in microtome chuck adapters and are resistant to the chemical action of most histological solvents. These cassettes are designed to accept standard metal lids (cat.# M481) and will keep specimens in complete safety during processing. The slanted writing surface accepts markings easily, permitting sample identification throughout all stages of embedding and long afterwards when in archives. They are available in 11 colors. Each case contains 3 dispenser boxes of 500 cassettes.

Cat. #	Color	Qty/Cs
M480-2	White	1500
M480-3	Pink	1500
M480-4	Green	1500
M480-5	Yellow	1500
M480-6	Blue	1500
M480-7	Peach	1500
M480-8	Tan	1500
M480-9	Gray	1500
M480-10	Lilac	1500
M480-11	Orange	1500
M480-12	Aqua	1500




M481 Metal Lid for M480 Cassette

Cat. #	Description	Qty/Pk
M481	Metal Process Cover	25



M480SL
For ThermoFisher Printers



M480T
For Leica and Sakura Printers

M480SL & M480T Cassettes in QuickLoad™ Sleeves and Stacks


Made of acetal

The sleeved cassettes are especially made to be used with ThermoFisher printers. Cassettes with tape are to be used with Leica and Sakura Ink Jet printers. Molded from a special high density acetal polymer, they keep specimens safely submerged and are resistant to the chemical action of most solvents used in histology laboratories. The efficient flow-through slots maximize fluid exchange and ensure proper reagent drainage.

In Sleeves			In Stacks		
Cat. #	Color	Qty/Cs	Cat. #	Color	Qty/Cs
M480-2SL	White	750	M480-2T	White	2000
M480-3SL	Pink	750	M480-3T	Pink	2000
M480-4SL	Green	750	M480-4T	Green	2000
M480-5SL	Yellow	750	M480-5T	Yellow	2000
M480-6SL	Blue	750	M480-6T	Blue	2000
M480-7SL	Peach	750	M480-7T	Peach	2000
M480-8SL	Tan	750	M480-8T	Tan	2000
M480-9SL	Gray	750	M480-9T	Gray	2000
M480-10SL	Lilac	750	M480-10T	Lilac	2000
M480-11SL	Orange	750	M480-11T	Orange	2000
M480-12SL	Aqua	750	M480-12T	Aqua	2000

Fluo Pink, Fluo Green and Fluo Yellow are also available. Minimum quantities apply. Contact us for more details.






M490



M491

Anatomy of a HISTOSETTE® I




The snap off lid can be easily removed, slid forward and snapped into the locked position.

Because it has efficient flow-through slots, fluid exchange and proper drainage are maximized.

Larger unobstructed writing surface sloped at 30° angle.

Because it has a back mounted locking device, it will never open during processing.



M490

HISTOSETTE® I

Tissue Processing / Embedding Cassettes

Made of acetal

Disposable plastic cassettes hold tissue specimens during the embedding process, as well as in a storage file. Molded from a special high density polymer, these cassettes keep specimens safely submerged in liquid and are totally resistant to the chemical action of most histological solvents. The efficient flow-through slots maximize fluid exchange and ensure proper drainage. The one-piece integral lid eliminates the need for separate steel lids. Just snap apart and conveniently lock the lid into the base of the cassette. They can be opened or closed as required they can be open or closed and they always relock securely without danger of specimen loss. Anterior writing area at a 30° angle. Not suitable for automated printers.

Each case contains 3 dispenser boxes of 500 cassettes.



M491

HISTOSETTE® I

Biopsy Processing / Embedding Cassettes

Made of acetal

Disposable plastic cassettes similar to Series M490 but specially designed to hold biopsy specimens during the embedding process, as well as in a storage file. Anterior writing area at a 30° angle. Not suitable for automated printers.

Each case contains 3 dispenser boxes of 500 cassettes.

Tissue		
Cat. #	Color	Qty/Cs
M490-2	White	1500
M490-3	Pink	1500
M490-4	Green	1500
M490-5	Yellow	1500
M490-6	Blue	1500
M490-7	Peach	1500
M490-8	Tan	1500
M490-9	Gray	1500
M490-10	Lilac	1500
M490-11	Orange	1500
M490-12	Aqua	1500

Biopsy		
Cat. #	Color	Qty/Cs
M491-2	White	1500
M491-3	Pink	1500
M491-4	Green	1500
M491-5	Yellow	1500
M491-6	Blue	1500
M491-7	Peach	1500
M491-8	Tan	1500
M491-9	Gray	1500
M491-10	Lilac	1500
M491-11	Orange	1500
M491-12	Aqua	1500

Fluo Pink, Fluo Green and Fluo Yellow are also available. Minimum quantities apply. Contact us for more details.



M498 HISTOSETTE® I

Tissue Processing / Embedding Cassettes

Made of acetal

These cassettes are identical to Series M490 but the anterior writing area has a 45 ° instead of a 30 ° angle. The acute angle makes these cassettes more suitable to be used with some models of cassette labeling instruments.

These disposable plastic cassettes hold tissue specimens very efficiently during the embedding process, as well as in a storage file. Molded from acetal, these cassettes keep specimens safely submerged in liquid and are resistant to the chemical action of most histological solvents. The efficient flow-through slots maximize fluid exchange and ensure proper drainage. The one-piece integral lid eliminates the need for separate steel lids. Just snap apart and conveniently lock the lid into the base of the cassette. They can be opened or closed ***as required they can be open or closed and they always relock securely without danger of specimen loss.*** Each case contains 3 dispenser boxes of 500 cassettes.

M499 HISTOSETTE® I

Biopsy Processing / Embedding Cassettes

Made of acetal

This model of cassette is similar to Series M498 and is specially designed to hold biopsy specimens during the embedding process.

Each case contains 3 dispenser boxes of 500 cassettes.



Fluo Pink, Fluo Green and Fluo Yellow are also available.
Minimum quantities apply. Contact us for more details.

Tissue		
Cat. #	Color	Qty/Cs
M498-2	White	1500
M498-3	Pink	1500
M498-4	Green	1500
M498-5	Yellow	1500
M498-6	Blue	1500
M498-7	Peach	1500
M498-8	Tan	1500
M498-9	Gray	1500
M498-10	Lilac	1500
M498-11	Orange	1500
M498-12	Aqua	1500

Biopsy		
Cat. #	Color	Qty/Cs
M499-2	White	1500
M499-3	Pink	1500
M499-4	Green	1500
M499-5	Yellow	1500
M499-6	Blue	1500
M499-7	Peach	1500
M499-8	Tan	1500
M499-9	Gray	1500
M499-10	Lilac	1500
M499-11	Orange	1500
M499-12	Aqua	1500

How to use a HISTOSETTE® I



Separate cover by folding it forward and backward.



Identify in front or sides.




Insert sample.



Slide cover in place and lock it.



After processing, remove cover by pulling on back tab.



M492 HISTOSETTE® II

Tissue Processing / Embedding Cassettes

Made of acetal

Disposable plastic cassettes hold tissue specimens during processing and embedding, as well as in storage. Molded from a special high-density polymer, these patented cassettes keep specimens safely submerged in liquid and are resistant to the chemical action of most histological solvents. **The efficient flow-through slots maximize fluid exchange and ensure proper drainage.** The one-piece integral lid eliminates the need for separate steel lids. They can be opened and closed ***as often as necessary*** and will always relock securely without danger of specimen loss.

The anterior writing area is slanted at a 45 ° angle to make the cassette more suitable to be used with certain types of cassette labeling instruments. Available in 11 colors. Each case contains 3 dispenser boxes of 500 cassettes.



Compatible with
most cassette printers

Tissue		
Cat. #	Color	Qty/Cs
M492-2	White	1500
M492-3	Pink	1500
M492-4	Green	1500
M492-5	Yellow	1500
M492-6	Blue	1500
M492-7	Peach	1500
M492-8	Tan	1500
M492-9	Gray	1500
M492-10	Lilac	1500
M492-11	Orange	1500
M492-12	Aqua	1500

Biopsy		
Cat. #	Color	Qty/Cs
M493-2	White	1500
M493-3	Pink	1500
M493-4	Green	1500
M493-5	Yellow	1500
M493-6	Blue	1500
M493-7	Peach	1500
M493-8	Tan	1500
M493-9	Gray	1500
M493-10	Lilac	1500
M493-11	Orange	1500
M493-12	Aqua	1500


Fluo Pink, Fluo Green and Fluo Yellow are also available.
Minimum quantities apply. Contact us for more details.

M493 HISTOSETTE® II

Biopsy Processing / Embedding Cassettes

Made of acetal



Similar to M492 but specially designed to hold biopsy specimens during the embedding process, as well as in a storage cabinet. Anterior writing area is at a 45 ° angle to make the cassette more suitable to be used with certain types of cassette labeling instruments. Each case contains 3 dispenser boxes of 500 cassettes.



How to use a HISTOSETTE® II

Unique front hinge design allows single handed cassette manipulation.

Unobstructed front writing surface sloped at 45 ° makes the cassette compatible with most automated cassette labeling machines.



Back mounted locking device securely holds the lid in place and will never open during processing.

Efficient flow-through slots in lid and base maximize fluid exchange, thereby ensuring proper drainage.



M485

HISTOSETTE® II

Tissue Processing / Embedding Cassettes

Made of acetal

Most convenient for cassette labeling instruments since covers are already removed from cassettes and are packaged separately in the case. Disposable plastic cassettes hold tissue specimens during processing and embedding, as well as in storage. Molded from a special high-density polymer, these patented cassettes keep specimens safely submerged in liquid and are resistant to the chemical action of most histological solvents. **The efficient flow-through slots maximize fluid exchange and ensure proper drainage.**

The anterior writing area is slanted at a 45 ° angle to make the cassette more suitable to be used with certain types of cassette labeling instruments. Available in 11 colors. Each case contains 2 dispenser boxes of 500 cassettes and 1 dispenser box of 1000 covers.

Tissue		
Cat. #	Color	Qty/Cs
M485-2	White	1000
M485-3	Pink	1000
M485-4	Green	1000
M485-5	Yellow	1000
M485-6	Blue	1000
M485-7	Peach	1000
M485-8	Tan	1000
M485-9	Gray	1000
M485-10	Lilac	1000
M485-11	Orange	1000
M485-12	Aqua	1000

Biopsy		
Cat. #	Color	Qty/Cs
M486-2	White	1000
M486-3	Pink	1000
M486-4	Green	1000
M486-5	Yellow	1000
M486-6	Blue	1000
M486-7	Peach	1000
M486-8	Tan	1000
M486-9	Gray	1000
M486-10	Lilac	1000
M486-11	Orange	1000
M486-12	Aqua	1000




Cassettes and lids
packaged separately



Compatible with
all cassette printers

Fluo Pink, Fluo Green and Fluo Yellow are also available.
Minimum quantities apply. Contact us for more details.



M486

HISTOSETTE® II

Biopsy Processing / Embedding Cassettes

Made of acetal


Most convenient for cassette labeling instruments but specially designed to hold biopsy specimens during the embedding process, as well as in a storage file. The covers are already removed from cassettes and are packaged separately in the case.

Anterior writing area is at a 45 ° angle to make the cassette more suitable to be used with certain types of cassette labeling instruments. Each case contains 2 dispenser boxes of 500 cassettes and 1 dispenser box of 1000 covers.


How to use a HISTOSETTE® II




Identify sample on
either the front
labeling area.




Press on back tab to
open cover.




Insert sample.



Close cover and press
front of lid to lock it in place.



To open, lift tab at the
back of cassette with
forefinger while depressing
the center of lid with thumb.



To remove the cover,
pivot the lid forward
and it will disconnect
automatically.

**M502****MICROSETTE® I****Biopsy Processing / Embedding Cassettes**

Made of acetal

Biopsy pads are no longer necessary with these innovative disposable plastic biopsy cassettes with a large compartment measuring 25 x 30 mm. Perfect even with needle biopsies since mesh holes have a diameter of only 0.26 mm while still allowing for maximum fluid exchange and drainage. Molded from a special high density polymer, these patented cassettes keep specimens safely submerged in liquid and are resistant to most histological solvents. They can be opened and closed as often as necessary and they always relock securely without danger of specimen loss. Available in 11 colors. Anterior writing area is at a 45° angle to make the cassette more suitable to be used with certain types of cassette labeling instruments.


Each case contains four dispenser boxes of 250 cassettes with covers assembled.




Compatible with
most cassette printers

Anatomy of a MICROSETTE™ I

Back mounted locking device securely holds the lid in place and will never open during processing.



No biopsy pads necessary.



Cover and base have over 2000, 0.26 mm square openings to maximize fluid exchange and ensure proper drainage. Excellent fluid exchange through slots.

Air vents allowing more efficient filling with paraffin.

Biopsy Cassette with one compartment

Cat. #	Color	Qty/Cs
M502-2	White	1000
M502-3	Pink	1000
M502-4	Green	1000
M502-5	Yellow	1000
M502-6	Blue	1000
M502-7	Peach	1000
M502-8	Tan	1000
M502-9	Gray	1000
M502-10	Lilac	1000
M502-11	Orange	1000
M502-12	Aqua	1000

Biopsy Cassette with six compartments

Cat. #	Color	Qty/Cs
M503-2	White	1000
M503-3	Pink	1000
M503-4	Green	1000
M503-5	Yellow	1000
M503-6	Blue	1000
M503-7	Peach	1000
M503-8	Tan	1000
M503-9	Gray	1000
M503-10	Lilac	1000
M503-11	Orange	1000
M503-12	Aqua	1000

Fluo Pink, Fluo Green and Fluo Yellow are also available. Minimum quantities apply. Contact us for more details.


M503**MICROSETTE® I****Biopsy Processing / Embedding Cassettes**

Made of acetal

This model can hold up to six tissue specimens, each one placed in its own 7 x 12 mm (1/4 x 7/16 in.) compartment, numbered from 1 to 6. No biopsy pads necessary. Cover and base have over 2000 square openings to maximize fluid exchange and ensure proper drainage. Approximately 170 holes (each having a diameter of 0.26 mm) per compartment.

Each case contains four dispenser boxes of 250 cassettes with covers assembled.





M505 & M506 UNISETTE® I

Biopsy & Tissue Processing / Embedding Cassettes

Made of acetal

Molded from a special high density polymer, these cassettes keep specimens safely submerged in liquid and are resistant to the chemical action of most histological solvents. The efficient flow-through slots maximize fluid exchange and ensure proper drainage. The one-piece integral lid eliminates the need for separate steel lids. The snap-latch and hinge-lock design of the UNISETTE™ prevent early separation of base and lid and allow for one-hand operation. Lids can be opened and closed **as required they can be open or closed and they always relock securely without danger of specimen loss**. Anterior writing area is at a 35° angle.

M506 Series is similar to M505 but specially designed to hold biopsy specimens during the embedding process.

Available in 11 colors. Each case contains 3 dispenser boxes of 500 cassettes.

Tissue			Biopsy		
Cat. #	Color	Qty/Cs	Cat. #	Color	Qty/Cs
M505-2	White	1500	M506-2	White	1500
M505-3	Pink	1500	M506-3	Pink	1500
M505-4	Green	1500	M506-4	Green	1500
M505-5	Yellow	1500	M506-5	Yellow	1500
M505-6	Blue	1500	M506-6	Blue	1500
M505-7	Peach	1500	M506-7	Peach	1500
M505-8	Tan	1500	M506-8	Tan	1500
M505-9	Gray	1500	M506-9	Gray	1500
M505-10	Lilac	1500	M506-10	Lilac	1500
M505-11	Orange	1500	M506-11	Orange	1500
M505-12	Aqua	1500	M506-12	Aqua	1500



Compatible with
most cassette printers

Have you ever considered
The SLIMSETTE™ ?

See M509



Fluo Pink, Fluo Green and Fluo Yellow are also available.
Minimum quantities apply. Contact us for more details.

How to use a UNISETTE™



The cover of the Unisetite is attached to the back of the base in an open position.



To close, simply tilt cover forward, holding it by the front tab.



Push cover down until a click is heard.



To open, tilt cover backward. It will remain attached to the base and can be closed again.



M507 MICROMESH™

NO BIOPSY PADS REQUIRED

Biopsy Processing / Embedding Cassettes

Made of acetal

This version of the Micromesh™ offers 1676 square openings (0.38 mm) allowing for a greatly improved fluid exchange and drainage without having to use biopsy pads. Large anterior slots in both cassette and cover ensure that the Micromesh™ biopsy cassette will sink rapidly. A large square compartment measuring 27 mm is perfect even for needle biopsies. The cover does not protrude above the cassette, a great space saving feature allowing more cassettes to be stacked in automatic labeling machines and tissue processors.

Molded from a special high density polymer, these patented cassettes keep specimens safely submerged in liquid and are resistant to the chemical action of most histological solvents. They can be opened and closed as required they can be open or closed and they always relock securely without danger of specimen loss. Anterior writing area is at a 45° angle to make the cassette more suitable to be used with automated cassette printers. Available in 11 colors.

Each case contains four dispenser boxes of 250 cassettes with covers assembled.

FOR
IVD
USE



Compatible with
most cassette printers


Biopsy Cassette with one compartment

Cat. #	Color	Qty/Cs
M507-2	White	1000
M507-3	Pink	1000
M507-4	Green	1000
M507-5	Yellow	1000
M507-6	Blue	1000
M507-7	Peach	1000
M507-8	Tan	1000
M507-9	Gray	1000
M507-10	Lilac	1000
M507-11	Orange	1000
M507-12	Aqua	1000

Biopsy Cassette with four compartments

Cat. #	Color	Qty/Cs
M508-2	White	1000
M508-3	Pink	1000
M508-4	Green	1000
M508-5	Yellow	1000
M508-6	Blue	1000
M508-7	Peach	1000
M508-8	Tan	1000
M508-9	Gray	1000
M508-10	Lilac	1000
M508-11	Orange	1000
M508-12	Aqua	1000

Anatomy of MICROMESH™ Cassettes



A recessed cover is a great space saving feature, allowing more cassettes to be stacked in automatic labeling machines and storage drawers

Fluo Pink, Fluo Green and Fluo Yellow are also available.
Minimum quantities apply. Contact us for more details.

M508 MICROMESH™

Biopsy Processing / Embedding Cassettes

Made of acetal

This model is similar to Series M507 but cassettes have four square compartments each measuring 13 mm. Cover and base have about 1676 square openings maximizing fluid exchange and ensuring proper drainage.

Each case contains four dispenser boxes of 250 cassettes with covers assembled.

With four compartments
each measuring 13 mm





M509 SLIMSETTE™

Tissue Processing / Embedding Cassettes

Made of acetal

The SLIMSETTE™, a new generation emerging in the Simport® Histology Family. More compact, easier to use and more efficient than ever. Similar to the design of the M507 model, it incorporates a unique recessed cover, a great space saving feature allowing more cassettes to be stacked in automatic labeling machines and in storage cabinets.

Molded from a special high density polymer, these cassettes keep specimens safely submerged in liquid and are resistant to the chemical action of most histological solvents. The SLIMSETTE™ ensures efficient fluid exchange and drainage thanks to 114 openings each measuring 1 x 5 mm. They can be opened and closed as required they can be open or closed and they always relock securely without danger of specimen loss. Large labeling areas for easy identification. The anterior writing area is slanted at a 45 ° angle. Available in 11 colors.

Each case contains 3 dispenser boxes of 500 cassettes with covers assembled.

Dimensions: 41 x 28.5 x 6 mm H (1 5/8 x 1 1/8 x 1/4 in. H)



Tissue Cassette with one compartment

Cat. #	Color	Qty/Cs
M509-2	White	1500
M509-3	Pink	1500
M509-4	Green	1500
M509-5	Yellow	1500
M509-6	Blue	1500
M509-7	Peach	1500
M509-8	Tan	1500
M509-9	Gray	1500
M509-10	Lilac	1500
M509-11	Orange	1500
M509-12	Aqua	1500

Biopsy Cassette with one compartment


Cat. #	Color	Qty/Cs
M510-2	White	1500
M510-3	Pink	1500
M510-4	Green	1500
M510-5	Yellow	1500
M510-6	Blue	1500
M510-7	Peach	1500
M510-8	Tan	1500
M510-9	Gray	1500
M510-10	Lilac	1500
M510-11	Orange	1500
M510-12	Aqua	1500

Tissue Cassette with four compartments

Cat. #	Color	Qty/Cs
M511-2	White	1500
M511-3	Pink	1500
M511-4	Green	1500
M511-5	Yellow	1500
M511-6	Blue	1500
M511-7	Peach	1500
M511-8	Tan	1500
M511-9	Gray	1500
M511-10	Lilac	1500
M511-11	Orange	1500
M511-12	Aqua	1500



Compatible with all cassette printers



M510 SLIMSETTE™

Biopsy Processing / Embedding Cassettes

Made of acetal

Similar to M509 but specially designed to hold biopsy specimens during the embedding process, as well as in a storage cabinet. Anterior writing area is at a 45 ° angle to make the cassette more suitable to be used with certain types of cassette labeling instruments. This biopsy model ensures efficient fluid exchange and drainage thanks to 392 openings



M511 SLIMSETTE™

Tissue Processing / Embedding Cassettes

Made of acetal

This model is similar to Series M509 but cassettes have four square compartments each measuring 13 mm. Cover and base have openings maximizing fluid exchange and ensuring proper drainage. Each case contains 3 dispenser boxes of 500 cassettes with covers assembled.

**M515****SWINGSETTE™****Tissue Processing / Embedding Cassettes**

Made of acetal

These disposable plastic cassettes hold tissue specimens during processing and embedding, as well as in storage. Molded from a special high-density polymer, they keep specimens safely submerged in liquid and are resistant to the chemical action of most histological solvents. The efficient flow-through slots maximize fluid exchange and ensure proper drainage.

This new model differs by the special hinge holding the base and cover together. This hinge allows cassettes to be opened and closed *as often as necessary*. The cover can be removed and re-inserted easily without danger of specimen loss.

Available in 11 colors.

Each case contains 3 dispenser boxes of 500 cassettes.



**A great alternative to
the Richard-Allan and
Surgipath Cassettes**

Anatomy of a SWINGSETTE™

Large tab for
convenient and easy
opening of lid

Secure locking
device holds the lid
in place during
processing

Unique hinge allows
opening and closing of
pre-installed lid as
often as needed

45° angle
writing area

Tissue		
Cat. #	Color	Qty/Cs
M515-2	White	1500
M515-3	Pink	1500
M515-4	Green	1500
M515-5	Yellow	1500
M515-6	Blue	1500
M515-7	Peach	1500
M515-8	Tan	1500
M515-9	Gray	1500
M515-10	Lilac	1500
M515-11	Orange	1500
M515-12	Aqua	1500

Biopsy		
Cat. #	Color	Qty/Cs
M516-2	White	1500
M516-3	Pink	1500
M516-4	Green	1500
M516-5	Yellow	1500
M516-6	Blue	1500
M516-7	Peach	1500
M516-8	Tan	1500
M516-9	Gray	1500
M516-10	Lilac	1500
M516-11	Orange	1500
M516-12	Aqua	1500

Fluo Pink, Fluo Green and Fluo Yellow are also available.
Minimum quantities apply. Contact us for more details.

M516**SWINGSETTE™****Biopsy Processing / Embedding Cassettes**

Made of acetal

These biopsy cassettes are similar to series M515 but especially designed to hold biopsy specimens during the processing / embedding process as well as in storage cabinets.

Each case contains three dispenser boxes of 500 cassettes.





M517 SWINGSETTE™

Tissue Processing / Embedding Cassettes

Made of acetal

Most convenient for cassette labeling instruments since the covers are packaged separately in the case. These disposable plastic cassettes hold tissue specimens during processing and embedding, as well as in storage. Molded from a special high-density polymer, they keep specimens safely submerged in liquid and are resistant to the chemical action of most histological solvents. The efficient flow-through slots maximize fluid exchange and ensure proper drainage.

This model differs by the special hinge that holds the base and cover together. This hinge allows cassettes to be opened and closed *as often as necessary*. The cover can be removed and re-inserted easily without danger of specimen loss. Available in 11 colors.

Each case contains two dispenser boxes of 500 cassettes and one dispenser box of 1000 covers.

Tissue			Biopsy		
Cat. #	Color	Qty/Cs	Cat. #	Color	Qty/Cs
M517-2	White	1000	M518-2	White	1000
M517-3	Pink	1000	M518-3	Pink	1000
M517-4	Green	1000	M518-4	Green	1000
M517-5	Yellow	1000	M518-5	Yellow	1000
M517-6	Blue	1000	M518-6	Blue	1000
M517-7	Peach	1000	M518-7	Peach	1000
M517-8	Tan	1000	M518-8	Tan	1000
M517-9	Gray	1000	M518-9	Gray	1000
M517-10	Lilac	1000	M518-10	Lilac	1000
M517-11	Orange	1000	M518-11	Orange	1000
M517-12	Aqua	1000	M518-12	Aqua	1000

Fluo Pink, Fluo Green and Fluo Yellow are also available.
Minimum quantities apply. Contact us for more details.




Cassettes and lids
packaged separately



Compatible with
all cassette printers

FOR
IVD
USE



M518 SWINGSETTE™

Biopsy Processing / Embedding Cassettes

Made of acetal

Most convenient for cassette labeling instruments since the covers are packaged separately in the case. These biopsy cassettes are similar to the M517 Series but are especially designed to hold biopsy specimens during the processing / embedding process as well as in storage cabinets.

Each case contains two dispenser boxes of 500 cassettes and one box of 1000 covers.


Close-up of the hinge



M515-M516 Series
The cover is connected
to the base.



M517-M518 Series
The special hinge makes it easy to attach lid to base.





M485SL & M486SL HISTOSETTE® II

Cassettes in QuickLoad™ Sleeves

Made of acetal


Suitable for hoppers accepting plastic sleeves, these cassettes will load in cassette labeling instruments in one simple operation. Save time and money with these convenient sleeves of 75 unique Simport® cassettes. No more manual insertion, one cassette at a time. Just load the cassette sleeve in the hopper and you are ready for printing.

These specially designed cassettes belong to the world-wide known Simport® HISTOSETTE® II Series. Molded from a special high density polymer, these cassettes keep specimens safely submerged in liquid and are resistant to the chemical action of most histological solvents. The efficient flow-through slots maximize fluid exchange and ensure proper drainage.

Two types of cassettes are offered in order to suit your particular needs: a regular tissue cassette and a biopsy model designed to hold small biopsy samples securely during the embedding process. Anterior printing area is at a 45 ° angle, offering an unobstructed view of the writing surface and making the cassette perfectly suitable to be used with cassette labeling instruments.

The one-piece disposable plastic cover eliminates the need for reusable steel lids. It can be opened and closed *as often as necessary* and it always relocks without danger of specimen loss. Available in 11 colors.

Each case contains 10 sleeves and 750 covers.



Tissue

Cat. No	Color	Qty/Cs
M485-2SL	Blanc	750
M485-3SL	Rose	750
M485-4SL	Vert	750
M485-5SL	Jaune	750
M485-6SL	Bleu	750
M485-7SL	Pêche	750
M485-8SL	Tan	750
M485-9SL	Gris	750
M485-10SL	Lilas	750
M485-11SL	Orange	750
M485-12SL	Aqua	750

Biopsy

Cat. No	Color	Qty/Cs
M486-2SL	Blanc	750
M486-3SL	Rose	750
M486-4SL	Vert	750
M486-5SL	Jaune	750
M486-6SL	Bleu	750
M486-7SL	Pêche	750
M486-8SL	Tan	750
M486-9SL	Gris	750
M486-10SL	Lilas	750
M486-11SL	Orange	750
M486-12SL	Aqua	750

Fluo Pink, Fluo Green and Fluo Yellow are also available.

Minimum quantities apply. Contact us for more details.



M482 & M483 HISTOSETTE® II


Cassettes in E-Z Load™ Stacks

Made of acetal

This stack of cassettes allows you to load the Shur/Mark® Cassette Labeling Instruments in one simple operation. Save time and money with these convenient stacks of 50 unique Simport® cassettes. No more manual insertion, one cassette at a time. Just load the stack in the hopper, cut and remove the holding tie and you are ready for printing.

These specially designed cassettes belong to the world-wide known Simport® Histosette® II design. Molded from a special high density polymer, these cassettes keep specimens safely submerged in liquid and are totally resistant to the chemical action of most histological solvents. The efficient flow-through slots maximize fluid exchange and ensure proper drainage.

Two types of cassettes are offered in order to suit your particular needs: a regular tissue cassette and a biopsy model designed to hold small biopsy samples securely during the embedding process. Anterior printing area is at a 45 ° angle.



TISSUE CASSETTES (In stacks of 50)

Base Cat. #	Cover Cat. #	Color	Qty/Cs
M482-2SL	M482-2LI	Blanc	2000
M482-3SL	M482-3LI	Rose	2000
M482-4SL	M482-4LI	Vert	2000
M482-5SL	M482-5LI	Jaune	2000
M482-6SL	M482-6LI	Bleu	2000
M482-7SL	M482-7LI	Pêche	2000
M482-8SL	M482-8LI	Tan	2000
M482-9SL	M482-9LI	Gris	2000
M482-10SL	M482-10LI	Lilas	2000
M482-11SL	M482-11LI	Orange	2000
M482-12SL	M482-12LI	Aqua	2000

BIOPSY CASSETTES (In stacks of 50)

Base Cat. #	Cover Cat. #	Color	Qty/Cs
M483-2SL	M483-2LI	Blanc	2000
M483-3SL	M483-3LI	Rose	2000
M483-4SL	M483-4LI	Vert	2000
M483-5SL	M483-5LI	Jaune	2000
M483-6SL	M483-6LI	Bleu	2000
M483-7SL	M483-7LI	Pêche	2000
M483-8SL	M483-8LI	Tan	2000
M483-9SL	M483-9LI	Gris	2000
M483-10SL	M483-10LI	Lilas	2000
M483-11SL	M483-11LI	Orange	2000
M483-12SL	M483-12LI	Aqua	2000

Fluo Pink, Fluo Green and Fluo Yellow are also available. Minimum quantities apply. Contact us for more details.



These cassettes can be used in this automatic labeling instrument.

Cassettes bases and lids are ordered separately



Cassettes for Primera® Printer

Made of acetal


Disposable plastic cassettes hold tissue specimens during processing and embedding, as well as in storage. Molded from a special high-density polymer, these cassettes keep specimens safely submerged in liquid and are resistant to the chemical action of most histological solvents. The efficient flow-through slots maximize fluid exchange and ensure proper drainage. As required they can be opened or closed and they always relock securely without danger of specimen loss. Available in 11 colors.

Specially designed for the automated and manual Primera Cassette printers.




Tissue

Biopsy



QuickLoad™ Stacks for Robotic feed printer

HISTOSETTE® II




TISSUE

BIOPSY

45° Angle

UNISETTE™ II




TISSUE

BIOPSY

35° Angle












SLIMSETTE™ II

























TISSUE


BIOPSY

45° Angle

Automated feed printer		Manual feed printer		
50 stacks per case Case 2000		Individual Cassette 2 dispenser boxes of 500 cassettes and 1 dispenser box of 1000 covers / Case 1000		
TISSUE	BIOPSY	TISSUE	BIOPSY	Color
Product #	Product #	Product #	Product #	
M385-2T	M386-2T	M385-2	M386-2	 White
M385-3T	M386-3T	M385-3	M386-3	 Pink
M385-4T	M386-4T	M385-4	M386-4	 Green
M385-5T	M386-5T	M385-5	M386-5	 Yellow
M385-6T	M386-6T	M385-6	M386-6	 Blue
M385-7T	M386-7T	M385-7	M386-7	 Peach
M385-8T	M386-8T	M385-8	M386-8	 Tan
M385-9T	M386-9T	M385-9	M386-9	 Gray
M385-10T	M386-10T	M385-10	M386-10	 Lilac
M385-11T	M386-11T	M385-11	M386-11	 Orange
M385-12T	M386-12T	M385-12	M386-12	 Aqua

Automated feed printer		Manual feed printer		
25 stacks per case Case 1000		Individual Cassette 3 dispenser boxes of 500 cassettes / Case 1500		
TISSUE	BIOPSY	TISSUE	BIOPSY	Color
Product #	Product #	Product #	Product #	
M405-2T	M406-2T	M405-2	M406-2	 White
M405-3T	M406-3T	M405-3	M406-3	 Pink
M405-4T	M406-4T	M405-4	M406-4	 Green
M405-5T	M406-5T	M405-5	M406-5	 Yellow
M405-6T	M406-6T	M405-6	M406-6	 Blue
M405-7T	M406-7T	M405-7	M406-7	 Peach
M405-8T	M406-8T	M405-8	M406-8	 Tan
M405-9T	M406-9T	M405-9	M406-9	 Gray
M405-10T	M406-10T	M405-10	M406-10	 Lilac
M405-11T	M406-11T	M405-11	M406-11	 Orange
M405-12T	M406-12T	M405-12	M406-12	 Aqua

Automated feed printer		Manual feed printer		
50 stacks per case Case 2000		Individual Cassette 3 dispenser boxes of 500 cassettes / Case 1500		
TISSUE Product #	BIOPSY Product #	TISSUE Product #	BIOPSY Product #	Color
M409-2T	M410-2T	M409-2	M410-2	 White
M409-3T	M410-3T	M409-3	M410-3	 Pink
M409-4T	M410-4T	M409-4	M410-4	 Green
M409-5T	M410-5T	M409-5	M410-5	 Yellow
M409-6T	M410-6T	M409-6	M410-6	 Blue
M409-7T	M410-7T	M409-7	M410-7	 Peach
M409-8T	M410-8T	M409-8	M410-8	 Tan
M409-9T	M410-9T	M409-9	M410-9	 Gray
M409-10T	M410-10T	M409-10	M410-10	 Lilac
M409-11T	M410-11T	M409-11	M410-11	 Orange
M409-12T	M410-12T	M409-12	M410-12	 Aqua



QuickLoad™ Stacks
for Robotic feed printer

M407 MICROMESH™ for Primera® Printer

Made of acetal

The MICROMESH™ offers 1676 square openings (0.38 mm) allowing for a greatly improved fluid exchange without having to use biopsy pads. Large anterior and posterior slots in both cassette and cover ensure that the MICROMESH™ Biopsy Cassette will sink rapidly. A large square compartment with a side measuring 27 mm is perfect even for needle biopsies. The cover does not protrude above the cassette, a great space saving feature.

Molded from a special high density polymer, these cassettes keep specimens safely submerged in liquid and are resistant to the chemical action of most histological solvents. The MICROMESH™ mesh ensures efficient fluid exchange and drainage. The one-piece integral lid eliminates the need for separate steel lids.


They can be opened and closed and always relock securely without danger of specimen loss. Anterior writing area is at a 45° angle.



Individual Cassette
for manual feed printer

Anatomy of Micromesh™ Cassettes

A recessed cover is a great space saving feature, allowing more cassettes to be stacked in automatic labeling machines and storage drawers




M407/M407T MICROMESH™

Individual Cassette for
manual feed printer
4 dispenser boxes of 250 /
Case 1000

QuickLoad™ Stacks for
robotic feed printer
50 stacks per case /
case 2000

TISSUE Product #	BIOPSY Product #	Color
M407-2	M407-2T	White
M407-3	M407-3T	Pink
M407-4	M407-4T	Green
M407-5	M407-5T	Yellow
M407-6	M407-6T	Blue
M407-7	M407-7T	Peach
M407-8	M407-8T	Tan
M407-9	M407-9T	Gray
M407-10	M407-10T	Lilac
M407-11	M407-11T	Orange
M407-12	M407-12T	Aqua

Primera® is a registered trademark of Primera Technology Incorporated



M507SL MICROMESH™

Biopsy Cassettes in QuickLoad™ Sleeves

Made of acetal

Suitable for hoppers accepting plastic sleeves, these cassettes will load in cassette labeling instruments in one simple operation. Just load the 75-cassette QuickLoad™ sleeve in the hopper and you are ready for printing.

The Micromesh™ offers 0.38 mm square openings, and large anterior slots allowing for increased fluid exchange and faster sinking in liquids. No biopsy pads are necessary. These cassettes feature a recessed cover and a large square 27 mm compartment, perfect even for needle biopsies.

Molded from a special high density polymer, these cassettes keep specimens safely submerged in liquid and are resistant to the chemical action of histological solvents. The lid can be opened and closed *as often as necessary* and always relocks securely without danger of specimen loss. Anterior writing area is at a 45 ° angle to make the cassette more suitable to be used with automated cassette printers. Available in 11 colors. Each case contains 10 sleeves.

FOR
IVD
USE



Compatible with
ThermoFisher
cassette printers

Lids are pre-mounted
on cassettes

NO BIOPSY PADS
REQUIRED



Biopsy

Cat. #	Color	Qty/Cs
M507-2SL	White	750
M507-3SL	Pink	750
M507-4SL	Green	750
M507-5SL	Yellow	750
M507-6SL	Blue	750
M507-7SL	Peach	750
M507-8SL	Tan	750
M507-9SL	Gray	750
M507-10SL	Lilac	750
M507-11SL	Orange	750
M507-12SL	Aqua	750


Fluo Pink, Fluo Green and Fluo Yellow are also available. Minimum quantities apply. Contact us for more details.

Anatomy of a MICROMESH™

A recessed cover is a great space saving feature, allowing more cassettes to be stacked in automatic labeling machines and storage drawers



Base and cover together have 1676 square openings (0.38 mm) allowing for a greatly improved fluid exchange



M509

Lid already installed
on cassette base



Compatible with
ThermoFisher
cassette printers

M509SL SLIMSETTE™

Tissue Cassettes in QuickLoad™ Sleeves

Made of acetal

The transparent sleeve allows viewing of cassettes in order to confirm there is no jam in the sleeve during the printing process.

This model is similar to the M507 cassette but is intended for tissue processing/embedding procedures. The efficient flow-through slots maximize fluid exchange and ensure proper drainage. Lids can be opened and closed as required they can be open or closed and they always relock securely without danger of specimen loss. When ready for the embedding process, covers can conveniently be snapped off and discarded.

Anterior writing area at a 45 ° angle. Available in 11 colors. Each case contains 10 sleeves of 75 cassettes.



Tissue		
Cat. #	Color	Qty/Cs
M509-2SL	White	750
M509-3SL	Pink	750
M509-4SL	Green	750
M509-5SL	Yellow	750
M509-6SL	Blue	750
M509-7SL	Peach	750
M509-8SL	Tan	750
M509-9SL	Gray	750
M509-10SL	Lilac	750
M509-11SL	Orange	750
M509-12SL	Aqua	750

Biopsy		
Cat. #	Color	Qty/Cs
M510-2SL	White	750
M510-3SL	Pink	750
M510-4SL	Green	750
M510-5SL	Yellow	750
M510-6SL	Blue	750
M510-7SL	Peach	750
M510-8SL	Tan	750
M510-9SL	Gray	750
M510-10SL	Lilac	750
M510-11SL	Orange	750
M510-12SL	Aqua	750

Fluo Pink, Fluo Green and Fluo Yellow are also available.
Minimum quantities apply. Contact us for more details.

M510SL SLIMSETTE™


Biopsy Cassettes in QuickLoad™ Sleeves

Made of acetal

Suitable for hoppers accepting plastic sleeves, these biopsy cassettes will load in cassette labeling instruments in one simple operation. Just load the 75-cassette QuickLoad™ sleeve in the hopper and you are ready for printing.

Molded from a special high density polymer, these cassettes keep specimens safely submerged in liquid and are totally resistant to the chemical action of most histological solvents. The lid can be opened and closed as often as necessary and always relocks securely without danger of specimen loss. Anterior writing area is at a 45° angle to make the cassette more suitable to be used with automated cassette printers. Available in 11 colors.

Cassettes are packaged in sleeves of 75. Each case contains 10 sleeves.



M510



M517SL SWINGSETTE™

Tissue Cassettes in QuickLoad™ Sleeves

Made of acetal

Suitable for hoppers accepting plastic sleeves, these tissue cassettes will load in cassette labeling instruments in one simple operation. Just load the 75-cassette QuickLoad™ sleeve in the hopper and you are ready for printing.

Molded from a special high density polymer, these cassettes keep specimens safely submerged in liquid and are resistant to the chemical action of histological solvents. The lid can be opened and closed *as often as necessary* and always relocks securely without danger of specimen loss. Anterior writing area is at a 45° angle to make the cassette more suitable to be used with automated cassette printers. Available in 11 colors.

Each case contains 10 sleeves of 75 cassettes and 10 bags of 75 lids.



Tissue			Biopsy		
Cat. #	Color	Qty/Cs	Cat. #	Color	Qty/Cs
M517-2SL	White	750	M518-2SL	White	750
M517-3SL	Pink	750	M518-3SL	Pink	750
M517-4SL	Green	750	M518-4SL	Green	750
M517-5SL	Yellow	750	M518-5SL	Yellow	750
M517-6SL	Blue	750	M518-6SL	Blue	750
M517-7SL	Peach	750	M518-7SL	Peach	750
M517-8SL	Tan	750	M518-8SL	Tan	750
M517-9SL	Gray	750	M518-9SL	Gray	750
M517-10SL	Lilac	750	M518-10SL	Lilac	750
M517-11SL	Orange	750	M518-11SL	Orange	750
M517-12SL	Aqua	750	M518-12SL	Aqua	750

Fluo Pink, Fluo Green and Fluo Yellow are also available. Minimum quantities apply. Contact us for more details.



Compatible with
ThermoFisher
cassette printers

M518SL SWINGSETTE™


Biopsy Cassettes in QuickLoad™ Sleeves

Made of acetal

These disposable plastic cassettes hold biopsy specimens during processing and embedding, as well as in storage. Molded from a special high-density polymer, they keep specimens safely submerged in liquid and are resistant to the chemical action of histological solvents. The efficient flow-through slots maximize fluid exchange and ensure proper drainage.


This new model differs by the special hinge holding the base and cover together. This hinge allows cassettes to be opened and closed *as often as necessary*. The cover can be removed and re-inserted easily without danger of specimen loss. Available in 11 colors.


Each case contains 10 sleeves of 75 cassettes and 10 bags of 75 lids.



Have you ever considered
Our Paraffin Block Mailer?

See M477-6





M492T HISTOSETTE® II

Tissue Cassettes in QuickLoad™ Stacks

Made of acetal

Specially made to be used with Leica and Sakura Ink Jet printers. They keep specimens safely submerged and are resistant to the chemical action of most solvents used in histology laboratories. The efficient flow-through slots maximize fluid exchange and ensure proper reagent drainage. The one-piece disposable plastic cover is pre-installed on each cassette and eliminates the need for reusable steel lids. The cover can be opened and closed as often as necessary and will always relock, reducing the possibility of specimen loss.

Cassettes are packaged in stacks of 40. Each case contains 50 stacks for a total of 2000 cassettes. Choose from 11 standard colors. The anterior printing area is at a 45 ° angle and offers an unobstructed view of the writing surface.



Compatible with LEICA and SAKURA cassette labeling instruments

Cat. #	Color	Qty/Cs
M492-2T	White	2000
M492-3T	Pink	2000
M492-4T	Green	2000
M492-5T	Yellow	2000
M492-6T	Blue	2000
M492-7T	Peach	2000
M492-8T	Tan	2000
M492-9T	Gray	2000
M492-10T	Lilac	2000
M492-11T	Orange	2000
M492-12T	Aqua	2000



M493T HISTOSETTE® II

Biopsy Cassettes in QuickLoad™ Stacks

Made of acetal

Histosette II® Biopsy Cassettes are similar to M492 Tissue Cassettes but specially designed to hold biopsy specimens. They have an attached lid that opens from the back of the cassette. The closed lids can be opened many times, always relocking securely. The anterior printing area is at a 45 ° angle and offers an unobstructed view of the writing surface.

Cassettes are packaged in stacks of 40. Each case contains 50 stacks for a total of 2000 cassettes. Choose from 11 standard colors.

Cat. #	Color	Qty/Cs
M493-2T	White	2000
M493-3T	Pink	2000
M493-4T	Green	2000
M493-5T	Yellow	2000
M493-6T	Blue	2000
M493-7T	Peach	2000
M493-8T	Tan	2000
M493-9T	Gray	2000
M493-10T	Lilac	2000
M493-11T	Orange	2000
M493-12T	Aqua	2000



Cat. #	Color	Qty/Cs
M505-2T	White	1000
M505-3T	Pink	1000
M505-4T	Green	1000
M505-5T	Yellow	1000
M505-6T	Blue	1000
M505-7T	Peach	1000
M505-8T	Tan	1000
M505-9T	Gray	1000
M505-10T	Lilac	1000
M505-11T	Orange	1000
M505-12T	Aqua	1000

M505T UNISETTE™

Tissue Cassettes in QuickLoad™ Stacks

Made of acetal

Molded from a special high density acetal, these tissue cassettes keep specimens safely submerged in liquid and are resistant to the chemical action of most solvents used in histology laboratories. The efficient flow-through slots maximize fluid exchange and ensure proper drainage.

The one-piece snap-latch and hinge-lock design prevents early separation of base and lid and allows one-hand operation. Lids are attached in an open position for easy filling, but can be opened or closed *as often as necessary* and will always relock securely. The tab on the left front side makes opening easy. Lids are easily removed by pulling sideways. The anterior writing area is at a 35° angle. Cassettes are packaged in stacks of 40. Each case contains 25 stacks for a total of 1000 cassettes. Choose from 11 standard colors.

FOR
IVD
USE



Compatible with LEICA
and SAKURA cassette
labeling instruments

Cat. #	Color	Qty/Cs
M506-2T	White	1000
M506-3T	Pink	1000
M506-4T	Green	1000
M506-5T	Yellow	1000
M506-6T	Blue	1000
M506-7T	Peach	1000
M506-8T	Tan	1000
M506-9T	Gray	1000
M506-10T	Lilac	1000
M506-11T	Orange	1000
M506-12T	Aqua	1000


M506T UNISETTE™

Biopsy Cassettes in QuickLoad™ Stacks

Made of acetal

The UNISETTE™ Biopsy Cassettes, made of acetal, are specially designed to hold biopsy specimens. One millimeter openings maximize fluid exchange and ensure proper drainage. The anterior writing area is at a 35° angle. The lids are attached, but arrive open for easy filling. There is a tab for opening on the left front side of the cassette lid. Lids are easily removed and will always relock securely.

Cassettes are packaged in stacks of 40. Each case contains 25 stacks for a total of 1000 cassettes. Choose from 11 standard colors.



M507T MICROMESH™

Biopsy Cassettes in QuickLoad™ Stacks

Made of acetal

The MICROMESH™ offers 1676 square openings (0.38 mm) allowing for a greatly improved fluid exchange without having to use biopsy pads. Large anterior and posterior slots in both cassette and cover ensure that the MICROMESH™ Biopsy Cassette will sink rapidly. A large square compartment with a side measuring 27mm is perfect even for needle biopsies. The cover does not protrude above the cassette, a great space saving feature.

Molded from a special high density polymer, these patented cassettes keep specimens safely submerged in liquid and are resistant to the chemical action of most histological solvents. The MICROMESH™ mesh ensures efficient fluid exchange and drainage. The one-piece integral lid eliminates the need for separate steel lids. They can be opened and closed ***as required they can be open or closed and they always relock securely without danger of specimen loss***. Anterior writing area is at a 45° angle to make the cassette more suitable to be used with automated cassette printers. Choose from 11 standard colors. Cassettes are packaged in stacks of 40. Each case contains 50 stacks for a total of 2000 cassettes.

Cat. #	Color	Qty/Cs
M507-2T	White	2000
M507-3T	Pink	2000
M507-4T	Green	2000
M507-5T	Yellow	2000
M507-6T	Blue	2000
M507-7T	Peach	2000
M507-8T	Tan	2000
M507-9T	Gray	2000
M507-10T	Lilac	2000
M507-11T	Orange	2000
M507-12T	Aqua	2000

FOR
IVD
USE



Compatible with LEICA
and SAKURA cassette
labeling instruments



M509T SLIMSETTE™

Tissue Cassettes in QuickLoad™ Stacks

Made of acetal

To be used with Leica and Sakura Ink Jet printers. Molded from a special high density acetal, they keep specimens safely submerged and are resistant to the chemical action of most solvents used in histology laboratories. The efficient flow-through slots maximize fluid exchange and ensure proper reagent drainage. The cover does not protrude above the cassette, a great space saving feature.

The anterior printing area is at a 45 ° angle and offers an unobstructed view of the writing surface.

The one-piece disposable plastic cover is pre-installed on each cassette and eliminates the need for reusable steel lids. The cover can be opened and closed ***as often as necessary*** and will always relock, reducing the possibility of specimen loss. Cassettes are packaged in stacks of 40. Each case contains 50 stacks for a total of 2000 cassettes. Choose from 11 standard colors.

Cat. #	Color	Qty/Cs
M509-2T	White	2000
M509-3T	Pink	2000
M509-4T	Green	2000
M509-5T	Yellow	2000
M509-6T	Blue	2000
M509-7T	Peach	2000
M509-8T	Tan	2000
M509-9T	Gray	2000
M509-10T	Lilac	2000
M509-11T	Orange	2000
M509-12T	Aqua	2000

**M510T****SLIMSETTE™****Biopsy Cassettes in QuickLoad™ Stacks**

Made of acetal

Slimsette™ Biopsy Cassettes are similar to M509 Tissue Cassettes but are specially designed to hold biopsy specimens. They have an attached lid that opens from the back of the cassette. The lids arrive closed but can be opened many times, always relocking securely. The anterior printing area is at a 45 ° angle and offers an unobstructed view of the writing surface.

Cassettes are packaged in stacks of 40. Each case contains 50 stacks for a total of 2000 cassettes. Choose from 11 standard colors.

Lid already installed on cassette base

Cat. #	Color	Qty/Cs
M510-2T	White	2000
M510-3T	Pink	2000
M510-4T	Green	2000
M510-5T	Yellow	2000
M510-6T	Blue	2000
M510-7T	Peach	2000
M510-8T	Tan	2000
M510-9T	Gray	2000
M510-10T	Lilac	2000
M510-11T	Orange	2000
M510-12T	Aqua	2000



Compatible with LEICA and SAKURA cassette labeling instruments

M517T & M518T
SWINGSETTE™**Tissue & Biopsy Cassettes in QuickLoad™ Stacks**


Made of acetal

Our latest model suited for the Leica and Sakura labellers. These cassettes will also load in those cassette labeling instruments in one simple operation. Save time and money with these convenient stacks of 40 Simport® cassettes. Just load the cassette sleeve in the magazine and you are ready for printing.

These specially designed cassettes belong to the world-wide known Simport® Swingsette™ design. They differ by the special hinge that holds the base and cover together. This hinge allows the Swingsette™ to be opened and closed *as often as necessary*.

Two types of cassettes are offered in order to suit your particular needs: a regular tissue cassette and a biopsy model designed to hold small biopsy samples securely during the embedding process. Anterior printing area is at a 45 ° angle, offering an unobstructed view of the writing surface and making the cassette perfectly suitable to be used with cassette labeling instruments. Choose from 11 standard colors.

TISSUE		BIOPSY	
Cat. #	Cat. #	Color	Qty/Cs
M517-2T	M518-2T	White	2000
M517-3T	M518-3T	Pink	2000
M517-4T	M518-4T	Green	2000
M517-5T	M518-5T	Yellow	2000
M517-6T	M518-6T	Blue	2000
M517-7T	M518-7T	Peach	2000
M517-8T	M518-8T	Tan	2000
M517-9T	M518-9T	Gray	2000
M517-10T	M518-10T	Lilac	2000
M517-11T	M518-11T	Orange	2000
M517-12T	M518-12T	Aqua	2000

**M512 MACROSETTE®****Processing / Embedding Cassettes (with lid)**

Made of acetal

Disposable plastic cassettes designed to hold larger tissue specimens during the embedding process, as well as in a storage cabinet. Dimensions are exactly the same as the ones of a regular HISTOSETTE® cassette but the MACROSETTE® is twice as high (13 mm).

The one-piece integral lid eliminates the need for separate steel lids. It can be opened and closed *as often as necessary* and will always relock securely without danger of specimen loss. Large labeling areas are located on three sides of the cassettes for your convenience. Each case contains three dispenser boxes of 250 cassettes.

Dimensions: 40.1 x 28.5 x 13 mm H (1 9/16 x 1 1/8 x 1/2 in. H)

Cat. #	Color	Qty/Cs
M512	White	750

M475-10 Disposable Deep Base Mold

Made of PVC

Designed especially for M512 MACROSETTE®. Thanks to the specially formulated plastic material, it offers excellent thermal exchange. It has a smooth interior finish and rounded corners facilitating specimen removal.

Cat. #	Size (mm)	Vol	Qty/Cs
M475-10	37 x 24 x 10	5 ml	500

**M474 Base Molds**

Made of Stainless Steel

Sizes fit most cassettes. Superior thermal exchange. These molds are for all applications in specimen embedding with all styles of Embedding Rings and Cassettes. Manufactured from high-quality stainless steel for optimal thermal conductivity, the molds have a highly polished surface for easy paraffin block removal. Well corners are rounded for optimal paraffin ribboning.

Cat. #	Size (mm)	Qty/Pk
M474-1	7 x 7 x 5	12
M474-2	15 x 15 x 5	12
M474-3	24 x 24 x 5	12
M474-4	30 x 24 x 5	12
M474-5	37 x 24 x 5	12

**M475 Disposable Base Molds**

Made of PVC

Simport® disposable base molds offer ease and convenience of operation. They are inexpensive enough to be discarded after use, yet strong enough to be reused. Thanks to the specially formulated plastic material, they offer excellent thermal exchange. They have a smooth interior finish and rounded corners facilitating specimen removal. Also, they are available in the same variety of sizes as metal molds and can be used with the same styles or types of cassettes and embedding rings. Each case contains two dispenser boxes of 500 base molds.

Cat. #	Size (mm)	Qty/Cs
M475-1	7 x 7 x 5	1000
M475-2	15 x 15 x 5	1000
M475-3	24 x 24 x 5	1000
M475-4	30 x 24 x 5	1000
M475-5	37 x 24 x 5	1000